1L12

LED-Mikroskoplicht  
LED Microscope Light

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Art.-Nr. / Art.No.  
100-004825

modulLED®

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Vielen Dank

This high quality product has been designed and developed to require a minimum of maintenance. Please follow the guidelines below and read this manual in its entirety to ensure safe operation of the unit.

## Warning

- Please disconnect plug from the main socket before cleaning the device.
- Storing and using the device only in dry rooms.

## Safety Information

- To avoid eye damage, don't look directly into the light source when turned on.

## Symbols

The following symbols are used in this manual.

**CAUTION!** This symbol indicates a potential risk and alerts you to proceed with caution.

**Important:** The evaluation of eye safety occurs according to the standard CIE S009/E:2002 („photobiological safety of lamps and lamp systems“). Within the risk grouping system of this CIE standard, the LEDs specified in this data sheet fall into the „low risk“ group (relating to devices in the visible spectrum with an exposure time of 100s). Under real circumstances (for exposure time, eye pupils, observation distance), it can be assumed that no endangerment to the eye exists from these devices. As a matter of principle, however, it should be mentioned that intense light sources have a high secondary exposure potential due to their blinding effect. As is also true when viewing other bright light sources (e.g. headlights), temporary reduction in visual acuity and afterimages can occur, leading to irritation, annoyance, visual impairment, and even accidents, depending on the situation.

## Intended Use

Professional LED IL1-Series are especially designed for the application in research and industry. Modern LED technology makes this products to a professional and efficient alternative to conventional halogen illuminating systems. The 3W Power Led's have an average life of circa 50.000 operating hours. The color temperature is with typical 5.500K similar to natural daylight and therefore ideal for umpteen professional application.

The “modul-LED” concept from StarLight allows a flexible and shoot exchange of the LED Moduls or the optional application from different optics and colors. The light intensity can be infinitely variable regulated process-controlled with an analog potentiometer. The products from IL-Serie are optional according to the IP67 watertight and are suitable for extrem industrial insets.

### IL12

IL12 two-armed

### Art.-No.

100-004825

## Package Contents

### Description

IL12 additional the moduls  
Power Supply Plug, 24V with industry plug  
User Manual

### Quantity

1  
1  
1

## Getting Started

The IL12 is a plug & play device very easy to use.

Please check before use whether ...

...the power cord is connected.

...the brightness is set to the desired value.

## Description of Buttons and Switches

### a) Electronic Dimming

For the continuous dimming of the light and for the Complete switch on / off of the ring light.

### b) Operation Control


The blue LED glows as soon as the ring light has been switched on.

### c) LED Moduls

A flexible and shoot exchange of the LED Moduls or the optional application from different optics and colors

### d) Connection from the power adaptor

Please connect correctly before use.

## Troubleshooting

The information in the table below can help to resolve the possible operating problems quickly.

Problem	Causes	Solution
Power supply: No function at the device. No blue LED.	No network connection.	Please check, whether the power cord and all other cables are connected.

## Storage, Maintenance, Repair

The device is maintenance-free. Some points, however, you should take into consideration.

1. Please unplug the device before cleaning.
2. For the cleaning of the housing please use a smooth, semi-wet cloth, if necessary with a smooth household cleaning fluid.
3. No moisture should entering into the interior of the housing.
4. To avoid damages during the storing and transport please retain the packaging.
5. Repairs are only allowed by trained and experienced personnel. In case of improper handling the warranty expires.

## Technical Service

### Web Resources:

Visit our StarLight website at [www.starlight-sl.de](http://www.starlight-sl.de) for:

- Contact details of technical customer service.
- StarLight product overview and price list as well as infos about accessories and spare parts.
- additional product information and special offers.

### Contact:

For information or technical assistance contact your local distributor or directly StarLight Sensorik Lichtleitsysteme GmbH & Co. KG.

### Technical Specifications

Weight:	0,75kg
Diameter optics:	19,7 mm
Operating Voltage:	100-240 VAC / 50-60 Hz
LED colour:	Pur-white, neutral-white, warm-white, green, blue, red and amber-yellow
Colour temperature:	Depending on selection
LED life span:	about 50.000h (according to manufacturer information)
Light control:	continuous electronic dimming from 0%-100%
Material:	Aluminum anodized/ Plastic
Colour of housing:	black

### Accessories and Spare Parts for IL-Series

 	various LED modules	See price list
		Power adaptor plug 24V with industry plug - Art.No.: 100-005333 Operating voltage: 100-240V, 50-60Hz Output voltage: 24V DC

## Warranty

StarLight Sensorik Lichtleitsysteme GmbH & Co. KG warrants that this product will be free from defects in material and workmanship for a period of one (1) year from date of purchase. If a defect is present, StarLight will, at its option, repair, replace, or refund the purchase price of this product at no charge to you, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear.

For your protection, items being returned must be insured against possible damage or loss. This warranty shall be limited to the replacement of defective products.

It is expressly agreed that this warranty will be in lieu of all warranties of fitness and in lieu of the warranty of merchantability.

## Equipment Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste. Instead, it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you

## Notizen/Notes

ryf ag

Ryf AG  
Bettlachstrasse 2  
2540 Grenchen  
tel 032 654 21 00  
fax 032 654 21 09

[www.ryfag.ch](http://www.ryfag.ch)